

Prevantics[®] Device Swab

3.15% Chlorhexidine Gluconate and 70% Isopropyl Alcohol Swab intended to disinfect needleless access sites prior to use

Scrub Time

5
Seconds

Dry Time

5
Seconds

WHERE TO USE:

- Needleless access sites prior to use

REMEMBER:

- Designed for disinfection of needleless access sites between each use.
- Do not blot or wipe dry.
- Implementation support provided by PDI's dedicated sales and clinical team.



PDI

Two Nice-Pak Park, Orangeburg, NY 10962-1376 + 1-800-999-6423
+ Phone: (845) 365-1700 + Fax: (845) 398-5347 + pdihc.com

©2014 Professional Disposables International, Inc. All rights reserved.
PDI10147831